



Non-Invasive Halo

Care and Use Instructions

1. This orthosis is designed to immobilize the head and neck region of the spine.
2. This orthosis is used because a traditional halo could not be used, for example soft bones would not hold the pins in the skull.
3. This orthosis works by having “tacky” silicone material in the head band and chin pieces to hold the head. These areas need to be cleaned with soap and water frequently to prevent skin irritation.
4. Cleaning is possible by having an assistant slide a wash cloth between the chin strap and the person. Then make sure the person’s skin is dry. This cleaning technique works best with the patient lying flat on their back in bed.
5. This orthosis is usually worn for approximately 3 months, until the bone heals. Frequent follow-ups with the physician will keep you informed on how things are going and the particular details for your case. If there are problems with the brace please call MOPC at 800 598 6672.